

Baby Mosquito

► Item: 19-732-002 PK /5



Active tip
thickness: 0.3 mm

Regular Mosquito
is 0.5 mm

40%
thinner
THAN REGULAR
MOSQUITO BUR

FOR USE
WITH
WATER
ONLY

USE THIS BUR:

- 1 With invisalign to ARS tight and small areas
- 2 To relief minor crowding in non-extraction cases
- 3 To eliminate black triangles
- 4 To enhance teeth contacts

► Shorter shank
to access
difficult to
reach areas

IOS

International Orthodontic
Services, PC

A BRAND
BUILT ON
TRUST

12811 Capricorn St. Stafford, TX 77477 USA
Toll Free USA and Canada 1888-10S-8882
P +1 832 342 9487 F +1 281 340 4764 E info@iosortho.com